



## Distinct Pregnancy Test Frequently Asked Questions (FAQ)

### Q: How does the test work?

**A:** The Distinct Pregnancy Test detects a hormone in your urine that your body produces during pregnancy (hCG-human chorionic gonadotropin). The amount of pregnancy hormone increases as pregnancy progresses.

### Q: How soon after I suspect that I am pregnant can I take the test?

**A:** You can test your urine as early as the first day of your missed period. You can perform the test anytime of the day; however, if you are pregnant, first morning urine contains the most pregnancy hormone.

### Q: Do I have to test with first morning urine?

**A:** Although you can test at any time of the day, your first morning urine is usually the most concentrated of the day and would have the most hCG in it.

### Q: How accurate is the test?

**A:** A clinical evaluation was conducted comparing the results obtained using the Distinct Pregnancy Test to another commercially available pregnancy test. The clinical trial included 100 urine specimens: both assays identified 50 positive and 50 negative results. The results demonstrated >99% overall accuracy of the Distinct Pregnancy Test when compared to the other pregnancy test.

### Q: How sensitive is the test?

**A:** The Distinct Pregnancy Test detects hCG in urine at a concentration of 25 mIU/mL or greater.

### Q: What should I do if the result shows that I am pregnant?

**A:** If your test says you are pregnant you should see your doctor who can advise you on what steps you should take next.

### Q: How do I know that the test was run properly?

**A:** The appearance of a colored line at the Control Line region tells you that you followed the test procedure properly and the proper amount of urine was absorbed.

### Q: What should I do if the result shows that I am not pregnant?

**A:** If your test says you are not pregnant, it means no hCG has been detected in your urine.

If you do not start your period within a week of its due date, repeat the test with a new pregnancy test. If you receive the same result after repeating the test and you still do not get your period, you should see your doctor.

### Q: How does this test compare to other tests on the market

**A:** The Distinct Pregnancy Test performance is comparable to the market leaders equivalent products.

MANUFACTURER/PRODUCT	FIRST RESPONSE RAPID RESULT PREGNANCY TEST	CLEARBLUE REASSURANCE PREGNANCY TEST	DISTINCT PREGNANCY TEST
Display	Line Test	Line Test	Line Test
Speed of Results	As Early as 1 Minute	As Early as 1 Minute	As Early as 1 Minute
When to Test	Day of Missed Period	Day of Missed Period	Day of Missed Period
Accuracy	99%	99%	99%



aconlabs.com

ACON Laboratories, Inc.  
10125 Mesa Rim Road, San Diego, CA 92121, U.S.A.  
Tel: 1.858.875.8000  
Fax: 1.858.200.0729  
Email: info@aconlabs.com