



# Early Detection Pregnancy Test

Detects Pregnancy up to 6 Days Early  
Quick Results

Highly Reliable and >99% Accurate

- *Wide Tip*
- *Ergonomic Handle*

## **AVAILABLE FOR OEM OPPORTUNITIES**

*Offering customers an affordable and reliable testing option!*



Highly Accurate



# Detection of Pregnancy

Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced by the developing placenta shortly after conception. hCG levels will increase throughout pregnancy. Specific hCG concentrations can be used to predict different stages of pregnancy.

## General Information

The Distinct Early Detection Pregnancy Test is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine. The test is sensitive enough to detect pregnancy up to 6 days early.

- Detection: Human chorionic gonadotropin
- Specimen: Urine

• Test Time: 1 - 3 Minutes

- Shelf life: 24 months
- Storage Temp.: 2 - 30 °C / 36 - 86°F

## Features




- Wide Absorbent Tip
- + or - Results Interpretation
- Protected Results Window
- Ergonomic Thumb Grip

## Test Procedure and Interpretation



Urine Stream



5~10s

OR



Collected Sample

### Testing in your Urine Stream

Hold the device by the Thumb Grip with the Absorbent Tip pointed downwards directly in the urine stream.


**Keep the tip in the stream for 5 -10 seconds.**

### Testing your Collected Urine Sample

Collect a sample of your urine in a clean, dry container. Hold the device by the Thumb Grip with the Absorbent Tip pointed downwards into the collected urine sample.

**Keep the tip submerged for 5-10 seconds.**

'Pregnant' result may be seen as early as 1 minute. **Wait the full 3 minutes to confirm a 'Not Pregnant' result.**



## Ordering Information



PRODUCT NAME	CATALOGUE	FORMAT	SPECIMEN	PACKAGE
Early Detection Pregnancy Test	L031-20192	Midstream	Urine	2 Tests/Kit; 3 Tests/Kit



aconlabs.com

ACON Laboratories, Inc.  
San Diego, CA 92121, U.S.A.  
Tel: 1.858.875.8000  
Fax: 1.858.200.0729  
Email: info@aconlabs.com

