

Deliver Value and Control



Self-monitoring of blood glucose is an important component of diabetes treatment.

But high costs put the benefits of frequent monitoring out of reach for many diabetics.

Answer the needs of your market with *On-Call® Plus*:

- Reliable
- Affordable
- Less Pain



Value and Control with *On-Call® Plus*



Starter kit contains 15 test strips
and all components needed to start testing

- Biosensor technology
- Only 1µL sample size
- Accurate results in just 10 seconds
- Flexibility to test on the forearm, palm or fingertip
- Advanced data management with 7, 14, 30-day averaging and PC download
- Affordability for frequent monitoring



For an interactive demonstration, please visit:
<http://www.aconlab.com/demo/G3%20flash/index.htm>

To order, please contact your ACON® distributor today.

Product Name	Catalog No.	Contents
<i>On-Call® Plus</i> Blood Glucose Monitoring System	G113-111 ✓	All components needed to start testing
<i>On-Call® Plus</i> Blood Glucose Meter	G113-211 ✓	Meter, carrying case, manual and control solution 1
	G113-214	Meter, lancing device, lancets, control solution 1, carrying case and manual
<i>On-Call® Plus</i> Blood Glucose Test Strips	G133-111 ✓	50 test strips (2 vials of 25)
	G133-114 ✓	30 test strips (2 vials of 15)
<i>On-Call® Plus</i> Glucose Control Solution	G123-311 ✓	One bottle of Control Solution 1 and one bottle of Control solution 2
<i>On-Call® Plus</i> Lancing Device	G123-211 ✓	Lancing device and insert
<i>On-Call® Plus</i> Lancets	G123-111 ✓	100 sterile lancets (4 bags of 25)

✓ CE Marked  0197

