

Mission™ Urinalysis Reagent Strips for Home Testing



With the *Mission™* Urinalysis Reagent Strips for Home Testing, you can test in the comfort of your own home!

- **Effective Testing and Monitoring**
- **Convenient**
- **Accurate**
- **Affordable**



Global Diagnostics for Local Markets™

Mission™ Urinalysis Reagent Strips for Home Testing provide accurate results in the comfort of your home

Effective Testing and Monitoring

- Test at home for possible health problems related to kidney and liver function, Urinary Tract Infection (UTI), diabetes, and other various diseases
- Help your doctor by monitoring frequently and detecting abnormalities before it becomes serious

Convenient


- Perform a simple urine test in the comfort of your home
- Test urine samples privately and as often as needed
- Easy to test and read
- Instant results after 2 minutes

Accurate


- Obtain reliable results by using the same strips that your doctor uses

Affordable


- Inexpensive testing without the need of laboratory instruments or home equipment



Step 1. Immerse strip into urine and immediately remove.



Step 2. Remove excess urine



Step 3. Compare the test pad(s) to the color blocks on the color chart.

Catalog No.	No. of Parameters	Test Name [√]	Test Purpose [∇]	Intended Use	Individually Pouched Strips per Kit	Strips per Canister [♠]	Parameters List													
							ASC	GLU	BIL	KET	SG	BLO	pH	PRO	URO	NIT	LEU			
U031-011	1	Glucose	Diabetes Screening and Monitoring	Prof/OTC	1, 3, 6	100		*												
		Ketone	Diabetes and/or Low Carb Diet Monitoring																	
		pH ^{NEW}	Monitor pH to help prevent Kidney stones and to help diagnose Vaginal Infection												*					
U031-021	2	Glucose/Ketone	Diabetes Monitoring	Prof/OTC	1, 3, 6	100		*		*										
		Nitrite/Leukocyte	Urinary Tract Infection (UTI) Screening												*	*				
		Bilirubin/Urobilinogen ^{NEW}	Liver Function Screening							*					*					
U031-031	3	Blood/Nitrite/Leukocyte ^{NEW}	Urinary Tract Infection (UTI) Screening	Prof/OTC	1, 3, 6	100						*				*	*			
U031-061	6	See Parameter List	General Health Checking	Prof/OTC	1, 3, 6	100			*	*	*	*	*	*	*	*	*	*	*	
U031-111	11	See Parameter List	General Health Checking and Physical Exam	Prof/OTC	1, 3, 6	100	*	*	*	*	*	*	*	*	*	*	*	*	*	

[√]Other OTC combinations available upon request. 3 new combinations available for Professional/OTC use on Jan 25, 2008.

[∇]For reference use only. The Test Purpose may depend on local medical practices.

[♠]Also available in canisters of 25, 50, and 150 strips.

CE Marked for sale in the European Community 

We also offer other rapid diagnostic and medical products for:

Lateral Flow Rapid Tests, Clinical Chemistry including Urinalysis, Immunoassay EIA/ELISA, Blood Glucose Monitoring Systems and more.

Contact us for worldwide distribution and custom manufacturing (OEM) opportunities

