

OptiLance™

Push-Button Safety Lancets

Post Acute Care & POC Testing



Product Name	OptiLance™	BD Microtainer®	ACCU-CHEK® Safe-T-Pro	ActiLance®	ProLance™	Unistik® 3	Assure Lance®	SurgiLance™
Marketed by	ACON	BD	Roche	HTL		Owen Mumford	Arkray	MediPurpose
Double Spring Action	✓	✓	✓	✓	✓	x	x	x
Spring Quality	Metal	Metal	Metal	Metal	Metal	Metal	Plastic	Metal
Gentle Grip for Patient Comfort	✓	x	x	x	x	x	x	x
Ergonomic Design	✓	✓	✓	✓	✓	x	x	x
Permanent Needle Retraction	✓	✓	✓	✓	✓	x	✓	✓
Needle Protection System	✓	✓	✓	✓	✓	✓	✓	x
Protective Cap	✓	✓	✓	✓	✓	✓	✓	✓

ACON Laboratories does not guarantee the accuracy or completeness of the information and data, makes no representations of warranties, express or implied, and assumes no liability associated with use or reliance thereof. All trademarks are the property of their respective owners. Customer specifications excerpted from product samples, labeling and marketing materials that were current as of March 2016.

For informational purposes only.
 ©2016 ACON Laboratories, Inc. San Diego, CA 92121
 PN 2160085001 3/16

