



On Call® Express



On Call® Vivid



On Call® Vivid Pal



On Call® Plus

No coding required	✓	✓	✓	Code Chip
Enzyme	<b>Glucose oxidase</b>	Glucose oxidase	Glucose oxidase	Glucose oxidase
Sample volume	<b>0.4 µL</b>	0.8 µL	0.8 µL	1.0 µL
Test time	<b>4 seconds</b>	5 seconds	5 seconds	10 seconds
Alternate site testing	<b>Palm, forearm</b>	Palm, forearm	Palm, forearm	Palm, forearm
Audible fill detection	✓	✓	✓	✓
Redosing if underfilled	✓	✓	✓	X
Testing reminder alarms	5	5	5	X
Hypo / hyper alarms	✓	✓	X	X
Ketone test reminder	✓	✓	X	X
Meal markers	✓	✓	✓	X
Memory	<b>300 tests</b>	500 tests	500 tests	300 tests
Data averaging	<b>7, 14, 30, 60, 90 days; Pre-meal, Post-meal, Total</b>	7, 14, 30, 60, 90 days; Pre-meal, Post-meal, Total	7, 14, 30, 60 and 90 days; Pre-meal, Post-meal, Total	7, 14, 30 days
Test strip ejector	✓	✓	X	X

All trademarks are the property of their respective owners. Competitor specifications excerpted from user manuals and test strip package inserts that were current as of September 2014. ACON Laboratories does not guarantee the accuracy or completeness of the information and data, makes no representations of warranties, express or implied, and assumes no liability associated with use or reliance thereof. For informational purposes only.