








On-Call® Stealth

Lancing Device

More Comfort for
your Patients



**Bright Yellow Release Button
Indicates Ready to Use**


-  Includes clear cap for alternate site testing
 -  11 depth settings provide the user with ultimate control over penetration depth for maximum comfort
 -  For use with universal lancets
 -  Designed for easy gripping and smooth release
 -  Lancet ejector to avoid accidental needle sticks
 -  Durable and long lasting for repeated use
 -  Easy draw with light spring mechanism
- Back Detail*



9 out of 10 people found On Call Stealth® to be less painful than their current lancing device*

Lancing Device Pain Evaluation Study

	1 = No Pain	2 = Slightly Painful	3 = Somewhat Painful	4 = Painful	5 = Very Painful
Accu-Chek® FastClix	54	34	7	4	1
On Call® Stealth™	51	36	10	2	1
OneTouch Delica™	34	44	17	4	1
ReliOn® (CanAm)	34	42	16	7	1
Nipro TRUEdraw™	9	45	20	17	9



Combined Response (2 responses per Subject)

On Call® Stealth™ vs. Accu-Chek® FastClix

- 28 Responses ranked On Call® Stealth™ as less painful
- 26 Responses ranked Accu-Chek™ as less painful
- 46 Responses Ranked them the same

No Statistically Different?

On Call® Stealth™ vs. OneTouch Delica™

- 40 Responses ranked On Call® Stealth™ as less painful
- 20 Responses ranked OneTouch Delica™ as less painful
- 40 Responses Ranked them the same

Yes (p<0.05)

Statistically Different?

1. 50 blindfolded subjects x 2 responses per subject

2. Lancing depth for each device normalized to the minimum depth required for a 0.4 µL sample from each subject

*On Call Stealth® Consumer Pain Study, 2016

All trademarks are the property of their respective owners.

ACON Laboratories does not guarantee the accuracy or completeness of the information and data, makes no representations of warranties, express or implied, and assumes no liability associated with use or reliance thereof. For informational purposes only.