

On-Call[®] Pro

Blood Glucose Meter

Post Acute



On-Call[®] Pro



Assure[®] Platinum



Assure[®] Prism Multi



Medline Evencare[®] G2



Medline Evencare[®] G3



EPS Bio Technology EASYMax[®] L



Nipro TRUE METRIX[™] PRO



ONETOUCH ULTRA[®] 2



ACCU-CHEK[®] PERFORMA

	On-Call [®] Pro	Assure [®] Platinum	Assure [®] Prism Multi	Medline Evencare [®] G2	Medline Evencare [®] G3	EPS Bio Technology EASYMax [®] L	Nipro TRUE METRIX [™] PRO	ONETOUCH ULTRA [®] 2	ACCU-CHEK [®] PERFORMA
Enzyme	Glucose oxidase	Glucose oxidase	Glucose oxidase	Glucose oxidase	Glucose oxidase	Glucose oxidase	GDH-FAD	Glucose oxidase	Mutant GDH-PQQ
Sample Volume	0.4 µL	0.5 µL	0.5 µL	1.0µL	0.8 µL	0.6 µL	0.5 µL	1.0 µL	0.6 µL
Test Time	4 Sec	7 Sec	5 Sec	6 Sec	6 Sec	5 Sec	4 Sec	5 Sec	5 Sec
Strip Port Direction	Faces patient	Faces patient	Faces patient	Faces patient	Facing operator, not patient	Faces patient	Facing operator, not patient	Faces patient	Facing operator, not patient
Strip Ejector	✓	Yes	Yes	No	Yes	No	Yes	No	No
QC Lockout Option	✓	Yes	No	No	Yes	No	No	No	No
Alternate Site Testing	✓ Palm, forearm	None	Palm, forearm, calf, thigh	Palm, forearm	Palm, forearm	Palm, forearm	Forearm	Palm, forearm	Palm, forearm, upper arm, thigh, calf
Audible Fill Detection	✓	Yes	Yes	Yes	Yes	Yes	Yes	No	No

All trademarks are the property of their respective owners. Competitor specifications excerpted from user manuals and test strip package inserts that were current as of February 2016.

ACON Laboratories does not guarantee the accuracy or completeness of the information and data, makes no representations of warranties, express or implied, and assumes no liability associated with use or reliance thereof. For informational purposes only.

©2016 ACON Laboratories, Inc. San Diego, CA 92121

PN 2160058003 10/16