

# Foresight<sup>®</sup> Tumor Markers EIA Test Kits AFP and CEA



Foresight<sup>®</sup> Tumor Markers EIA Test Kits offer high quality and accurate AFP and CEA testing at a cost-effective price!

- *Accurate and Efficient*
- *Easy to Perform*
- *Early Diagnosis and Continuous Monitoring*





# Foresight® Tumor Markers EIA Test Kits

## Accurate and Efficient

- Advance screening and diagnosis
- Reliable performance for quantitative results
- Dependable and cost-effective
- Analytical sensitivity as low as 0.6 ng/mL for AFP and 0.5 ng/mL for CEA

## Easy to Perform

- Room temperature incubations for more convenient assay
- Color-coded and numbered labeling for easy access
- Breakaway 8-well strips in 96 microwell plate stored in resealable foil pouch

## Early Diagnosis and Continuous Monitoring

- All reagents are ready to use to eliminate initial preparation
- Assay time as short as 45 minutes for faster results
- Facilitates early patient treatment decisions
- Monitors treatment progress



Feature	AFP <sup>✓</sup>	CEA <sup>✓</sup>
Catalog No.	I231-2011	I231-2021
Tests per Kit*	96	96
Type of Test	Quantitative	Quantitative
Principle	Double Antibody Sandwich EIA	Double Antibody Sandwich EIA
Specimen	Serum/Plasma	Serum/Plasma
Specimen Volume	25 µL	25 µL
Conjugate Volume	100 µL	100 µL
Substrate A/B Volume	50 µL/50 µL	50 µL/50 µL
Stop Solution Volume	50 µL	50 µL
Total Incubation Time	45 min (30 min/15 min)	45 min (30 min/15 min)
Incubation Temperature	Room Temperature	Room Temperature
Wavelength	450/630-700 nm or 450 nm	450/630-700 nm or 450 nm
Calibrator Range	0-400 ng/mL	0-120 ng/mL
Analytical Sensitivity	0.6 ng/mL	0.38 ng/mL

\* 480 and 48 tests per kit available soon

✓ CE Marked for sale in the European Community 

Other EIA products available:

HIV 1/2/O, HIV 1/2/O Antigen/Antibody, HBsAg, HBsAb, HBeAg, HBeAb, HBcAb, HCV, HEV IgM<sup>✓</sup>, Syphilis<sup>✓</sup>, ToRCH, TSH<sup>✓</sup>,  
*H. pylori* Antigen<sup>✓</sup> and *H. pylori* IgG<sup>✓</sup>

Upcoming EIA products including T3, T4, *H. pylori* IgA and IgM, Free PSA, Total PSA, HTLV I/II and HEV IgG.

We also offer other rapid diagnostic and medical products for:

Clinical Chemistry including Urinalysis, Blood Glucose Monitoring and Others.

Contact us for worldwide distribution and customer manufacturing (OEM) opportunities



ACON Laboratories, Inc., 10125 Mesa Rim Road, San Diego, CA 92121, USA.  
 Tel: 1-858-875-8000 • Fax: 1-858-200-0729 • E-mail: [info@aconlabs.com](mailto:info@aconlabs.com)  
 Please visit our website for details: [www.aconlabs.com](http://www.aconlabs.com)

Authorized Distributor