

Foresight®

EIA Test Kits for Infectious Diseases, Fertility, Tumor Markers, and Thyroid

Reliable, High Quality
Testing at Competitive
Prices!

- Dependable and Complete
- Efficient Testing
- Convenient



ACON®

Global Diagnostics for Local Markets™

Foresight® EIA Test Kits for Infectious Diseases, Fertility, Tumor Markers, and Thyroid

Dependable and Complete

- High sensitivity and specificity
- Consistent quality with long shelf life
- Multiple EIA tests available

Efficient Testing

- Total incubation time as short as 45 minutes for most kits and 15 minutes for Rapid hCG kit
- Color-coded and numbered reagents for easy handling
- Ready to use reagents shorten initial preparation

Convenient

- All necessary reagents provided
- Can be performed on automated processors
- Breakaway 8-well strips in 96-microwell plate in resealable foil pouch



Group Name	Product Name	Catalog No.	Detection	Principle	Type of Test	Total Incubation Time
Infectious Diseases	HIV 1/2/O Antibody	I231-1011	HIV-1 (p24, gp41), HIV-2, (gp36) and Subtype O antibodies	Double antigen sandwich EIA	Qualitative	60 min (30/20/10 min)
	HIV 1/2/O Antigen/Antibody	I231-1281	HIV-1 p24 Ag, HIV-1 (p24, gp41), HIV-2, (gp36) and Subtype O antibodies	Double antigen/antibody sandwich EIA	Qualitative	120 min (60/30/30 min)
	HTLV I/II	I231-1291	HTLV I/II IgG, IgA, and IgM antibodies	Double antigen sandwich EIA	Qualitative	90 min (30/30/30 min)
	HBsAg	I231-1021	Ad and ay subtypes of Hepatitis B surface antigens	Double antibody sandwich EIA	Qualitative	Standard: 70 min (60/10 min) Enhanced: 150 min (120/30 min)
	HBsAb	I231-1051	Hepatitis B surface antibody	Double antigen sandwich EIA	Qualitative	45 min (30/15 min)
	HBsAb Quantitative	I231-1271	Hepatitis B surface antibody	Double antigen sandwich EIA	Quantitative	45 min (30/15 min)
	HBeAg	I231-1061	Hepatitis B envelope antigen	Double antibody sandwich EIA	Qualitative	45 min (30/15 min)
	HBeAb	I231-1071	Hepatitis B envelope antibody	Competitive EIA	Qualitative	45 min (30/15 min)
	HBcAb	I231-1081	Hepatitis B core antibody	Competitive EIA	Qualitative	45 min (30/15 min)
	HBcAb IgG	I231-1311	IgG antibodies to Hepatitis B core antigen	Indirect EIA	Qualitative	70 min (30/30/10 min)
	HCV	I231-1031	Various subtypes of HCV antibodies using recombinant antigens specific for core, NS3, NS4, and NS5	Indirect EIA	Qualitative	70 min (30/30/10 min)
	HEV IgM	I231-1211 ✓	Hepatitis E IgM antibodies	Immunocapture EIA	Qualitative	75 min (30/30/15 min)
	Syphilis	I231-1041 ✓	Syphilis IgG, IgM, and IgA using recombinant antigens specific for TP15, TP17, and TP47	Double antigen sandwich EIA	Qualitative	75 min (60/15 min)
	<i>H. pylori</i> Antigen	I231-1231 ✓	<i>H. pylori</i> antigen	Double antibody sandwich	Qualitative/Quantitative	70 min (60/10 min)
	<i>H. pylori</i> IgG	I231-1241 ✓	IgG antibodies to <i>H. pylori</i>	Indirect EIA	Qualitative/Quantitative	70 min (30/30/10 min)
	Toxoplasma IgG	I231-1091 ✓	IgG antibodies to <i>Toxoplasma gondii</i>	Indirect EIA	Qualitative/Quantitative	70 min (30/30/10 min)
	Toxoplasma IgM	I231-1101 ✓	IgM antibodies to <i>Toxoplasma gondii</i>	Immunocapture EIA	Qualitative	70 min (30/30/10 min)
	Rubella IgG	I231-1111 ✓	IgG antibodies to Rubella	Indirect EIA	Qualitative/Quantitative	70 min (30/30/10 min)
	Rubella IgM	I231-1121 ✓	IgM antibodies to Rubella	Immunocapture EIA	Qualitative	70 min (30/30/10 min)
	CMV IgG	I231-1131 ✓	IgG antibodies to CMV	Indirect EIA	Qualitative/Quantitative	70 min (30/30/10 min)
CMV IgM	I231-1141 ✓	IgM antibodies to CMV	Immunocapture EIA	Qualitative	70 min (30/30/10 min)	
HSV 1 IgG	I231-1151 ✓	IgG antibodies to HSV 1	Indirect EIA	Qualitative	70 min (30/30/10 min)	
HSV 1 IgM	I231-1161 ✓	IgM antibodies to HSV 1	Immunocapture EIA	Qualitative	70 min (30/30/10 min)	
HSV 2 IgG	I231-1171 ✓	IgG antibodies to HSV 2	Indirect EIA	Qualitative	70 min (30/30/10 min)	
HSV 2 IgM	I231-1181 ✓	IgM antibodies to HSV 2	Immunocapture EIA	Qualitative	70 min (30/30/10 min)	
HSV 1/2 IgG	I231-1191 ✓	IgG antibodies to HSV 1/2	Indirect EIA	Qualitative	70 min (30/30/10 min)	
HSV 1/2 IgM	I231-1201 ✓	IgM antibodies to HSV 1/2	Immunocapture EIA	Qualitative	70 min (30/30/10 min)	
Fertility	Rapid hCG	I231-4011 ✓	hCG antigen	Double antibody sandwich EIA	Qualitative/Semi-Quantitative	15 min (10/5 min)
	hCG	I231-4051 ✓	hCG antigen	Double antibody sandwich EIA	Quantitative	75 min (60/15 min)
	LH	I231-4021 ✓	LH antigen	Double antibody sandwich EIA	Quantitative	75 min (60/15 min)
	Prolactin-PRL	I231-4041 ✓	PRL antigen	Double antibody sandwich EIA	Quantitative	75 min (60/15 min)
	FSH	I231-4031 ✓	FSH antigen	Double antibody sandwich EIA	Quantitative	75 min (60/15 min)
Tumor Marker	AFP	I231-2011 ✓	Alpha-Fetoproteine (AFP)	Double antibody sandwich EIA	Quantitative	45 min (30/15 min)
	CEA	I231-2021 ✓	Carcinoembryonic antigen (CEA)	Double antibody sandwich EIA	Quantitative	45 min (30/15 min)
	Total PSA	I231-2031 ✓	Total Prostate-Specific Antigen (PSA)	Double antibody sandwich EIA	Quantitative	45 min (30/15 min)
Thyroid	TSH	I231-3011 ✓	Thyroid Stimulating Hormone (TSH)	Double antibody sandwich EIA	Quantitative	80 min (60/20 min)
	Total T4	I231-3021 ✓	Total Thyroxine (T4)	Competitive EIA	Quantitative	75 min (60/15 min)
	Free T4	I231-3031 ✓	Free Thyroxine (FT4)	Competitive EIA	Quantitative	75 min (60/15 min)
	Total T3	I231-3041 ✓	Total Triiodothyronine (T3)	Competitive EIA	Quantitative	75 min (60/15 min)
	Free T3	I231-3051 ✓	Free Triiodothyronine (FT3)	Competitive EIA	Quantitative	75 min (60/15 min)

✓ CE Marked for sale in the European Community

We also offer other rapid diagnostic and medical products for:

Blood Glucose Monitoring Systems, Clinical Chemistry including Urinalysis, Molecular Diagnostic Reagents, and much more

Contact us for worldwide distribution and custom manufacturing (OEM) opportunities.

