

Foresight®

Fertility EIA Test Kits



Foresight® Fertility EIA Test Kits offer high quality and broad Fertility assays at an affordable price!

- *Reliable*
- *Easy to Perform*
- *Convenient*



Global Diagnostics for Local Markets™

Foresight® Fertility EIA Test Kits

Reliable

- High sensitivity equivalent to market leaders
- Broad Fertility assays ideal for detection and monitoring of developmental and reproductive disorders
- Reliable results at competitive pricing

Easy to Perform

- Same procedure for all Fertility assays except Rapid hCG
- Room temperature incubations for added convenience
- No additional reagent preparation necessary

Convenient

- Standardized color-coded and numbered labeling for easy access
- Breakaway 8-well strips in 96 microwell plate in resealable foil pouch
- Suitable for automated or manual testing depending on laboratory needs
- Shelf life of 15 months



Feature	Rapid hCG ✓	hCG ✓	LH ✓	Prolactin - PRL ✓	FSH ✓
Catalog No.	I231-4011	I231-4051	I231-4021	I231-4041	I231-4031
Tests per Kit	96	96	96	96	96
Type of Test	Qualitative Semi-Quantitative	Quantitative	Quantitative	Quantitative	Quantitative
Principle	Sandwich EIA	Sandwich EIA	Sandwich EIA	Sandwich EIA	Sandwich EIA
Specimen	Serum/Urine	Serum	Serum	Serum	Serum
Specimen Volume	20 µL	20 µL	50 µL	50 µL	50 µL
Conjugate Volume	100 µL	100 µL	100 µL	100 µL	100 µL
Substrate A/B Volume	50 µL	50 µL	50 µL	50 µL	50 µL
Stop Solution Volume	50 µL	50 µL	50 µL	50 µL	50 µL
Total Incubation Time	15 min (10/5 min)	75 min (60/15 min)	75 min (60/15 min)	75 min (60/15 min)	75 min (60/15 min)
Incubation Temperature	Room temperature	Room temperature	Room temperature	Room temperature	Room temperature
Wavelength	450 or 450/630-700 nm	450 or 450/630-700 nm	450 or 450/630-700 nm	450 or 450/630-700 nm	450 or 450/630-700 nm
Calibrator Range	0-250 mIU/mL	0-1000 mIU/mL	0-200 mIU/mL	0-200 ng/mL	0-100 mIU/mL
Analytical Sensitivity	0.68 mIU/mL	1.0 mIU/mL	0.5 mIU/mL	1.0 ng/mL	0.6 mIU/mL
Shelf Life	15 months	15 months	15 months	15 months	15 months

✓ CE Marked for sale in the European Community 

Other EIA products available:

HIV, HBV, HCV, HEV, Syphilis, HTLV, H. pylori, ToRCH, AFP, CEA, PSA, TSH, T3, T4

We also offer other rapid diagnostic and medical products for:

Clinical Chemistry including Urinalysis, Blood Glucose Monitoring and Others.

Contact us for worldwide distribution and custom manufacturing (OEM) opportunities.



ACON Laboratories, Inc., 10125 Mesa Rim Road, San Diego, CA 92121, USA • Tel: 1-858-875-8000 • Fax: 1-858-200-0729 • E-mail: info@aconlabs.com
Please visit our website for details: www.aconlabs.com