

SARS-CoV-2 IgG/IgM Rapid Test

Provides an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection.*



20 Years Experience



Reliable Performance



Easy to Use




User Friendly



CE Marked

SARS-CoV-2 IgG/IgM Rapid Test



Coronavirus disease 2019 (COVID-19) is an infectious disease caused by severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2).

Symptoms: Common symptoms include fever, cough, and shortness of breath. Other symptoms may include muscle pain, sputum production, diarrhea, sore throat, loss of smell, and abdominal pain. While the majority of cases result in mild symptoms, some progress to pneumonia and multi-organ failure.


General Information

- Detection: IgM and IgG antibodies to SARS-CoV-2
- Sensitivity: 99.1% (95% CI: 95.1%-100%)
- Specificity: 98.3% (95% CI: 96.2%-99.5%)
- Accuracy: 98.5% (95% CI: 96.9%-99.5%)
- Specimen: Serum/Plasma/Whole Blood from Venous or Fingerstick
- Test time: 15-20 minutes
- Shelf life: 12 months
- Storage temperature: 2-30°C

Materials Provided

- Test cassettes
- Disposable specimen applicator
- Buffer
- Package insert

Test Procedure and Interpretation



Ordering Information

Product Name	Catalog No.	Format	Specimen	Package
ACON SARS-CoV-2 IgG/IgM Rapid Test	L031-11711 ✓	Cassette	S/P/WB	25 Tests/Kit

✓ CE marked

*Results from the SARS-CoV-2 IgG/IgM Rapid Test should not be used as the sole basis for diagnosis.



aconlabs.com

ACON Laboratories, Inc.
5850 Oberlin Drive, #340, San Diego, CA 92121, USA
Tel: 1.858.875.8000
Fax: 1.858.200.0729
Email: info@aconlabs.com