

# CentriVet™ Strips: Urine Reagent Strips for Animal Urinalysis



## Reliable Performance

- Use with CentriVet™ Urine Analyzer or read visually
- Analytical sensitivity comparable to market leaders
- High quality color chart ensures accurate visual reading

## Quick and Easy Testing

- Accurate results for all parameters within 60 seconds
- CentriVet™ Urine Analyzer provides printed report for record keeping

## Multiple Packaging Options

- Available in 50 and 100 strips per canister – offers cost savings and convenience for low and high volume testing
- 2 year shelf life for strips in unopened canisters
- 3 month shelf life for strips in opened canisters



## Specifications

Features	Specifications
Test Parameters	LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU
Storage Temperature	2 - 30°C
Unopened Strip Shelf Life	24 months from date of manufacture
Opened Strip Canister Use Life	3 months

## Ordering Information

Product Order No.	Product Description
U031-108-ABN04A	CentriVet™ Urine Reagent Test Strips, 50 strips/canister
U031-108-ABN04B	CentriVet™ Urine Reagent Test Strips, 100 strips/canister