



Influenza A/B Ag Rapid Test Frequently Asked Questions (FAQ)

Q: I see a blue line in the control region. What does this mean?

A: A blue line in the control region indicates that the test was performed correctly and that the amount of sample used was sufficient. A blue line in the test region indicates a positive result for influenza A or B.

Q: Can I use the test if I have a fever or other symptoms?

A: Yes, you can use the test if you have a fever or other symptoms. The test is designed to detect influenza A and B antigens, which are present in the body during an active infection.

Q: How long does the test take to run?

A: The test takes approximately 15 minutes to run. Results are visible in the test region within 10 minutes. A blue line in the test region indicates a positive result for influenza A or B.

Q: Can I use the test if I am pregnant or breastfeeding?

A: Yes, you can use the test if you are pregnant or breastfeeding. The test is safe and accurate for all ages and conditions.

Q: I see a blue line in the test region. What does this mean?

A: A blue line in the test region indicates a positive result for influenza A or B.

Q: I see a blue line in the test region and a blue line in the control region. What does this mean?

A: A blue line in the test region and a blue line in the control region indicate a positive result for influenza A or B. This is the expected result for a positive test.

Q: I see a blue line in the control region but no line in the test region. What does this mean?

A: A blue line in the control region but no line in the test region indicates a negative result for influenza A or B. This could mean that you do not have an active infection or that the test was not performed correctly.

Q: How do I know that the test was run properly?

A: The appearance of a colored line at the Control Line region tells you that you followed the test procedure properly and the proper amount of sample was used.



aconlabs.com

ACON Laboratories, Inc.
5850 Oberlin Drive, # 340
San Diego CA 92121
Tel: 1.858.875.8000
Fax: 1.858.200.0729
Email: info@aconlabs.com