


## HBV Hepatitis B Virus Quantitative PCR Test Kit

Hepatitis virus types B and C are the most common cause of liver cirrhosis, cancer and viral hepatitis-related deaths. *Promotor*® HBV Hepatitis B Virus Quantitative PCR Test Kit is an in vitro PCR test for the quantification of HBV DNA in human serum or plasma.



### Test Performance

HBV Hepatitis B Virus Quantitative PCR Test Kit	Specifications
Package Size	32 tests/kit
Technology	Real-time PCR
Sample Types	Human serum or plasma
Genotypes	HBV genotypes of A - H
Sensitivity	20 IU/mL
Linear Range	50 IU/mL - $5.0 \times 10^8$ IU/mL
Storage Conditions	-20 ± 5°C

### Order Information

Product Name	Catalog No.	Package Size
HBV Hepatitis B Virus Quantitative PCR Test Kit	P131-1511	32 Tests/Kit

### Workflow



\* Automated extraction system can be used. It is recommended to use Promotor® NES-32 Nucleic Acid Extraction System and Promotor® Nucleic Acid (DNA) Extraction Kit manufactured by ACON

IVD

We also offer other diagnostic and medical products including:

Blood Glucose and other meters, Lateral Flow tests, Urinalysis products, POCT products, Immunoassays and more.

Contact us for worldwide distribution and custom manufacturing (OEM) opportunities.



[www.aconlabs.com](http://www.aconlabs.com)

**ACON Laboratories, Inc.**  
 5850 Oberlin Drive, #340, San Diego, CA 92121, USA  
 Tel: 1.858.875.8000  
 Fax: 1.858.200.0729  
 Email: [info@aconlabs.com](mailto:info@aconlabs.com)